

## Advanced endoscopic linear cutter stapler



*"X" curved tip*

- **Endoscopic Linear Cutting Stapler Reload ( Stepped - Staple ) For Single-use**
- **Stepped-staple Technology Delivers Graduated Compression, Optimizing Interaction Between The Tissue and Stapler**



## Stapler Specification



Stapler code	Total length	OD
BZS- S	280mm	12.5mm
BZS- M	379mm	12.5mm
BZS- L	479mm	12.5mm

## Reload Specification



Type	Cartridge type	Color	Staple height (mm)
45 series Staple length	45T	Brown	2.0 ; 2.5 ; 3.0
	45P	Purple	3.0 ; 3.5 ; 4.0
	45B	Black	4.0 ; 4.5 ; 5.0
	45XT	Brown	2.0 ; 2.5 ; 3.0
60 series Staple length	45XP	Purple	3.0 ; 3.5 ; 4.0
	60T	Brown	2.0 ; 2.5 ; 3.0
	60P	Purple	3.0 ; 3.5 ; 4.0
	60B	Black	4.0 ; 4.5 ; 5.0
	60XT	Brown	2.0 ; 2.5 ; 3.0
	60XP	Purple	3.0 ; 3.5 ; 4.0

- Light Weight Design
- Precise Articulation
- Forward / Backward compatibility
- Improved Unloading Mechanism
- One-handed Grasping Mechanism

**"X" The special curved tip is designed for better vascular anastomosis in thoracic surgery.**

