



High end Laparoscopic instruments



Laparoscopic Monopolar

Laparoscopic Bipolar forceps

Universal Monopolar & Bipolar


Suction and Irrigation



Models: OD: 5mm
BG-CF-G-5*180mm
BG-CF-G-5*330mm
BG-CF-G-5*360mm
BG-CF-G-5*420mm

- Many kinds of head are available.
- 360°rotating head,easy to locate during surgery.

Single use shaft to prevent cross infection




Gallbladder grasping forceps



Short scissors



Straight dissectors



Intestinal forceps



Mini gallbladder grasping forceps



Curved scissors



Curved dissector





Tooth grasping forceps



Models: OD: 5mm
BG-CF-T-5*180mm
BG-CF-T-5*330mm
BG-CF-T-5*360mm
BG-CF-T-5*420mm

- Smart choice for kinds of laparoscopic surgery.
- 360°rotating head,easy to locate during surgery.


Single use tip to prevent cross infection



Gallbladder grasping forceps



Short scissors



Straight dissectors



Intestinal forceps



Mini gallbladder grasping forceps



Curved scissors



Curved dissector



Tooth grasping forceps

Universal Monopolar & Bipolar

Monopolar



Monopolar coagulation forceps

Models: OD: 5mm
BG-DH-5*180mm
BG-DH-5*330mm
BG-DH-5*360mm
BG-DH-5*420mm

Hand control



Foot control

Models: OD: 5mm
BG-HIVB-5*330mm
BG-HIVB-5*450mm
BG-FIVB-5*330mm
BG-FIVB-5*450mm

Advanced Monopolar coagulation hook with suction

- Long suction tube design for deep surgical procedures
- Abdominal Cavity, Thoracic Cavity & More
- Hand and Foot Control options that meet the current trends and requirements of today's surgeons
- 45 degree scratch and suction tube with functionality that ensure minimally invasive, less tissue damaged procedures
- Suction functionality with valve switch, smoke removal & minimized pneumoperitoneum during surgery



Universal Monopolar coagulation hook

Models: OD: 5mm
BG-IVB-5*330mm
BG-IVB-5*450mm

Bipolar




Bipolar coagulation forceps

Models: OD: 5mm
BG-SH-5*180mm
BG-SH-5*330mm
BG-SH-5*360mm
BG-SH-5*420mm

Laparoscopic Electrocoagulator

Suction & Irrigation




Specification:

OD: 5mm	Length (mm)
BG-ELD-5	160,280,330,380,420
BG-ELS-5	160,280,330,380,420
BG-EWD-5	160,280,330,380,420
BG-EWS-5	160,280,330,380,420
BG-ELW-5	160,280,330,380,420
BG-EWW-5	160,280,330,380,420

Suction & Irrigation kits

Specification:

OD: 5mm	Length (mm)
BG-LD-5	160,280,330,380,420
BG-LS-5	160,280,330,380,420
BG-WD-5	160,280,330,380,420
BG-WS-5	160,280,330,380,420
BG-LW-5	280,330,420
BG-WW-5	280,330,420



- Set electrocoagulation; Electro-excision; Suction and irrigation as a whole, maximizing surgical efficiency
- With 7mm suction tube, easy remove cysts and broken tissue
- Fast switching valves with obvious color distinction, Purple for "Irrigation"; White for "Suction"
- Effective smoke evacuation
- Single use to avoid cross infection



Accessory

+ Single use Retrieval bag




Nylon
Length: 10 / 12 / 15mm



TPU
Bag diameter: 60 / 80 / 100mm

Models:
BG-NZQ-10*230ml
BG-NZQ-12*310ml
BG-NZQ-18*800ml
BG-TZQ-60*150ml
BG-TZQ-80*150ml
BG-TZQ-100*300ml

+ Single use Veress Needle



Models:
BG-QF-12*120mm
BG-QF-15*150mm

+ Single use Trocar

- Detachable sleeves, two layers of seals design secures minimal air leakage & smooth in & outs.
- The inverted triangle stability thread design on sleeve provides greater fixation to the abdominal wall.
- Wide range: diameter(5mm/10mm/12mm/15mm), length (70mm/100mm/150mm), etc.



Models:
Blade
BG-GD-CCQ
- 5/10/12/15mm
Bladeless
BG-FL-CCQ
- 5/10/12/15mm

+ Single use Ease Filter

Length: 250+20cm



Models: OD: 7mm
BG-GL

+ Compatible cables




Disposable bipolar



Disposable Mono



Reusable bipolar



Reusable Mono

Models: Length: 3m
BG-YC-1
Disposable
BG-CF-1
Reusable