


## General Surgical staplers



### linear cutting stapler

Product name	Model		Cutting length	Staple height	Color
	Handle	Reload			
Linear Cutting Stapler	BG - LC60B	BG - LCR60B	60mm	3.8mm	Blue
	BG - LC60G	BG - LCR60G	60mm	4.8mm	Green
	BG - LC80B	BG - LCR80B	80mm	3.8mm	Blue
	BG - LC80G	BG - LCR80G	80mm	4.8mm	Green
	BG - LC100B	BG - LCR100B	100mm	3.8mm	Blue
	BG - LC100G	BG - LCR100G	100mm	4.8mm	Green




### linear stapler

Product name	Model		Cutting length	Staple height	Color
	Handle	Reload			
Linear Stapler	BG - LS30B	BG - LSR30B	30mm	3.8mm	Blue
	BG - LS30G	BG - LSR30G	30mm	4.8mm	Green
	BG - LS45B	BG - LSR45B	45mm	3.8mm	Blue
	BG - LS45G	BG - LSR45G	45mm	4.8mm	Green
	BG - LS60B	BG - LSR60B	60mm	3.8mm	Blue
	BG - LS60G	BG - LSR60G	60mm	4.8mm	Green



### Hemorrhoidal stapler

Product name	Model	OD	Staple height	Reload color
Hemorrhoidal Stapler	BG - HS32	32mm	4.0mm	White
	BG - HS34	34mm	4.0mm	Green



## **Circular stapler**

- **Smaller head diameter**  
Smaller head diameter of anvil provides easier control for insertion and removal during placement.
- **Shorter anvil height**  
Short anvil height design benefits the kinking anatomy of esophagus for better solution.
- **Bigger staple contacting surface**  
Bigger staple contacting surface will greatly guide staples into the staple pocket which enables leak-resistant and hemostatic anastomosis.

Product name	Model	OD	Staple height	Reload color
Circular Stapler	BG - CS23	23mm	4.0mm	Dark blue
	BG - CS25	25mm	4.0mm	Light blue
	BG - CS29	29mm	4.0mm	Green
	BG - CS32	32mm	4.0mm	White

