

## Universal endoscopic linear cutter stapler

*Suitable for various open & minimally surgeries to perform transection, excision or anastomosis of tissue.*



- 360° rotation, with 0°, 22°, 45° articulation in right and left directions, meet the requirements of different parts of the operation.



- 3 More in line with the use of the design, improve the efficiency of anastomosis.

### Stapler Specification

Product Name	Model	Total length	Shaft length	OD (mm)
Universal Cutting Stapler	BG - ECL	479mm	260mm	12.4
	BG - ECM	379mm	160mm	12.4
	BG - ECS	289mm	70mm	12.4



## Specification



### Single use blade

Ensure maximum sharpness, avoid risk of cross-infection, reduce the possibility of cancer metastasis.




### Extra-thin reinforced anvil

Provide optimal tissue grasp force, ensure better staple formation, avoid over compression.

## Reload Specification

Product Name	Model	Cutting length	Staple height	Reload color
Articulating Reload	BG - ECR6048A	60mm	4.8mm	Green
	BG - ECR6040A	60mm	4.0mm	Golden
	BG - ECR6035A	60mm	3.5mm	Blue
	BG - ECR6025A	60mm	2.5mm	White
	BG - ECR4548A	45mm	4.8mm	Green
	BG - ECR4535A	45mm	3.5mm	Blue
	BG - ECR4525A	45mm	2.5mm	White
	BG - ECR3035A	30mm	3.5mm	Blue
	BG - ECR3025A	30mm	2.5mm	White
Straight Reload	BG - ECR6048	60mm	4.8mm	Green
	BG - ECR6040	60mm	4.0mm	Golden
	BG - ECR6035	60mm	3.5mm	Blue
	BG - ECR6025	60mm	2.5mm	White
	BG - ECR4548	45mm	4.8mm	Green
	BG - ECR4535	45mm	3.5mm	Blue
	BG - ECR4525	45mm	2.5mm	White
	BG - ECR3035	30mm	3.5mm	Blue
	BG - ECR3025	30mm	2.5mm	White



### B type Staple formation

Better staple formation, avoid bleeding effectively, at the same time to ensure the safety.

